



Product Datasheet

CHUK antibody (orb135280)


Description	Rabbit polyclonal antibody to CHUK
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, IP, WB
Immunogen	A synthetic peptide of human IKK α (NP_001269.3).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	85KDa
Uniprot ID	O15111
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200ICC1:50-200, IF/ICC: 1:50-200, IP: 1:20-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunofluorescence analysis of A549 cell...



Immunohistochemistry of paraffin-embedde...