



## Product Datasheet

### CRK antibody (orb135312)

|                            |                                                             |
|----------------------------|-------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal to CRK                                    |
| <b>Reactivity</b>          | Human, Mouse                                                |
| <b>Conjugation</b>         | Unconjugated                                                |
| <b>Tested Applications</b> | WB                                                          |
| <b>Immunogen</b>           | Recombinant fusion protein of human CRK/CRKII(NP_058431.2). |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.          |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.                 |
| <b>Note</b>                | For research use only                                       |
| <b>Isotype</b>             | IgG                                                         |
| <b>Clonality</b>           | Polyclonal                                                  |
| <b>MW</b>                  | 37kDa                                                       |
| <b>Uniprot ID</b>          | <b>P46108</b>                                               |
| <b>Dilution Range</b>      | WB: 1:500-2000, IHC-P: 1:50-200                             |
| <b>Expiration Date</b>     | 12 months from date of receipt.                             |



Western blot analysis of extracts of var...



Immunohistochemical staining of human li...