

## Product Datasheet

### DR3 (TNFRSF25) Human Over-expression Lysate (orb1353244)

|                        |                                                                                                         |
|------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Transient overexpression lysate of tumor necrosis factor receptor superfamily, member 25 (TNFRSF25),... |
| <b>Target</b>          | TNFRSF25                                                                                                |
| <b>Note</b>            | For research use only                                                                                   |
| <b>Uniprot ID</b>      | <b>Q93038</b>                                                                                           |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                          |