

Product Datasheet

Stathmin antibody (orb13536)

Description

Stathmin antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, ICC, IF, IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human Stathmin. Please contact us for the exact immunogen sequence. The peptide is available as orb374769.

Target

Stathmin

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

17 kDa

Uniprot ID

[P54227](#), [P13668](#), [P16949](#)

NCBI

[001145454](#), [001138926](#)

Entrez


[3925](#)

Dilution Range


WB: 2 µg/ml, IF/ICC: 1:100, IHC-P: 1:100

Expiration Date


12 months from date of receipt.



WB analysis of rat liver tissue using an...



Immunohistochemical staining of mouse br...



Western blot analysis of rat kidney tiss...