



Product Datasheet

LEF1 antibody (orb13540)


Description	LEF1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from human LEF-1 (331-399/399 aa)
Target	LEF-1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	44 kDa
Uniprot ID	Q9UJU2
Dilution Range	WB=1:500-2000, Flow-Cyt=1ug/test
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Paraffin-embedded rat brain tissue fixed...



IHC-P of mouse tissue (LEF1 antibody at ...)