

## Product Datasheet

### Estrogen Receptor beta2 antibody (orb135475)

**Description** Rabbit polyclonal antibody to Estrogen Receptor beta2

**Reactivity** Human, Mouse, Rat

**Conjugation** Unconjugated

**Tested Applications** WB

**Immunogen** Recombinant fusion protein of human Estrogen Receptor beta(NP\_001201832.1).

**Preservatives** PBS with 0.05% proclin300,50% glycerol,pH7.3.

**Concentration** 1 mg/ml

**Storage** Store at -20°C. Avoid freeze / thaw cycles.

**Note** For research use only

**Isotype** IgG


**Clonality** Polyclonal

**MW** 59KDa


**Uniprot ID** [Q92731](https://www.uniprot.org/entry/Q92731)

**Dilution Range** WB: 1:500-2000, IHC-P: 1:50-200

**Expiration Date** 12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemical staining of mouse lu...