



Product Datasheet

LRIG2 antibody (orb13548)

Description	LRIG2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human LRIG2 (701-800/1065 aa)
Target	LRIG2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	117 kDa
Uniprot ID	O94898
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



IHC-P of mouse kidney tissue using LRIG2...



Immunohistochemical analysis of formalin...