


Product Datasheet

FGFR1 antibody (orb135512)

Description	Rabbit polyclonal antibody to FGFR1
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human FGFR1(NP_075598.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	6kDa/16kDa/33kDa/55-95kDa
Uniprot ID	P11362
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...