



Product Datasheet

CD11c antibody (orb13554)


Description	CD11c antibody
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human CD11c. Please contact us for the exact immunogen sequence. The peptide is available as orb12913.
Target	CD11c
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	127 kDa
Uniprot ID	P20702, Q9QXH4
Entrez	3687
Dilution Range	IHC-P: 1: 100-600
Expiration Date	12 months from date of receipt.



Western blot analysis of rat lymph node ...



IHC-P staining of human ovarian cancer t...



Immunohistochemical staining of paraffin...