



## Product Datasheet

### CFH antibody (orb135638)

|                            |                                                                       |
|----------------------------|-----------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal to CFH                                              |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                     |
| <b>Conjugation</b>         | Unconjugated                                                          |
| <b>Tested Applications</b> | ICC, IF, IHC, WB                                                      |
| <b>Immunogen</b>           | Recombinant fusion protein of human Complement factor H(NP_000177.2). |
| <b>Preservatives</b>       | PBS with 0.05% proclin300,50% glycerol,pH7.3.                         |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.                           |
| <b>Note</b>                | For research use only                                                 |
| <b>Isotype</b>             | IgG                                                                   |
| <b>Clonality</b>           | Polyclonal                                                            |
| <b>MW</b>                  | 150kDa                                                                |
| <b>Uniprot ID</b>          | <b>P08603</b>                                                         |
| <b>Dilution Range</b>      | WB: 1:500-2000, IHC-P: 1:50-200                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                       |



Immunohistochemical staining of human li...



Immunohistochemical staining of human li...