



## Product Datasheet

### HLA-A antibody (orb135651)


|                            |                                                   |
|----------------------------|---------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to HLA-A               |
| <b>Reactivity</b>          | Human                                             |
| <b>Conjugation</b>         | Unconjugated                                      |
| <b>Tested Applications</b> | IHC, WB                                           |
| <b>Immunogen</b>           | A synthetic peptide of human HLA-A(NP_002107.3).  |
| <b>Preservatives</b>       | PBS with 0.05% proclin300,50% glycerol,pH7.3.     |
| <b>Concentration</b>       | 1 mg/ml                                           |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.       |
| <b>Note</b>                | For research use only                             |
| <b>Isotype</b>             | IgG                                               |
| <b>Clonality</b>           | Polyclonal                                        |
| <b>MW</b>                  | 45KDa                                             |
| <b>Uniprot ID</b>          | <b>P04439, P30443</b>                             |
| <b>Dilution Range</b>      | WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                   |



IHC-P image of human neurilemmoma tissue...



IHC-P staining of mouse kidney tissue us...



Immunohistochemical analysis of formalin...