


Product Datasheet

HOXA11 antibody (orb135679)

Description	Rabbit polyclonal to HOXA11
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human HOXA11(NP_005514.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	30KDa
Uniprot ID	P31270
Dilution Range	WB: 1:200-500
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of HeLa...