



Product Datasheet

IL8 antibody (orb135736)

Description	Rabbit polyclonal to IL8
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide of human IL8(NP_000575.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	10KDa
Uniprot ID	P10145
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse brain tissue.



Western blot analysis of extracts of A-431 cells.