



Product Datasheet

IL13 antibody (orb135744)


Description	Rabbit polyclonal to IL13
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human IL13 (NP_002179.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	13KDa
Uniprot ID	P35225
Entrez	3596
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat brain tissue.



Immunohistochemical staining of human tissue.



Western blot analysis of extracts of recombinant protein.