


Product Datasheet

MMP9 antibody (orb13583)


Description

MMP9 antibody


Species/Host	Rabbit
Reactivity	Bovine, Canine, Equine, Guinea pig, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human MMP9. Please contact us for the exact immunogen sequence. The peptide is available as orb374791.
Target	MMP9
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	78 kDa
Uniprot ID	P14780 , P41245
NCBI	004994 , 21 , 004985
Entrez	4318
Dilution Range	WB: 2 µg/ml, IHC-P:1:200, IF/ICC: 1:200
Expiration Date	12 months from date of receipt.



IHC-P staining of mouse liver tissue



WB analysis of Mouse spleen (Lane 1), Ra...



Western blot analysis of human muscle ti...