



Product Datasheet

LRP5 antibody (orb135855)


Description	Rabbit polyclonal to LRP5
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human LRP5 (NP_002326.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	180-200KDa
Uniprot ID	O75197
Entrez	4041
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human placenta tissue.



Immunohistochemical staining of mouse lung tissue.



Western blot analysis of extracts of various tissues.