



Product Datasheet

MAPT antibody (orb135878)


Description	Rabbit polyclonal antibody to MAPT.
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide of human Tau(NP_058519.3).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45-75KDa
Uniprot ID	P10636
Dilution Range	WB: 1:500-2000, IHC-P: 1:100-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse li...



Western blot analysis of extracts of HT-...



Immunohistochemical staining of mouse ki...