



Product Datasheet

MAX antibody (orb135881)


Description	Rabbit polyclonal to MAX
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human MAX (NP_002373.3).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	19kDa
Uniprot ID	P61244
Entrez	4149
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse brain tissue.



Immunofluorescence analysis of U2OS cells.



Western blot analysis of extracts of MCF7 cells.