


Product Datasheet

MCAM antibody (orb135887)

Description	Rabbit polyclonal to MCAM
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human CD146/CD146/MCAM(NP_006491.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	120KDa
Uniprot ID	P43121
Dilution Range	WB: 1:1000-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of HeLa...