



## Product Datasheet

### MEF2C antibody (orb135898)

|                            |                                                           |
|----------------------------|-----------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to MEF2C                       |
| <b>Reactivity</b>          | Human, Mouse, Rat                                         |
| <b>Conjugation</b>         | Unconjugated                                              |
| <b>Tested Applications</b> | IHC, WB                                                   |
| <b>Immunogen</b>           | A synthetic peptide of human MEF2C(NP_002388.2).          |
| <b>Preservatives</b>       | PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.          |
| <b>Concentration</b>       | 1 mg/ml                                                   |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.               |
| <b>Note</b>                | For research use only                                     |
| <b>Isotype</b>             | IgG                                                       |
| <b>Clonality</b>           | Polyclonal                                                |
| <b>MW</b>                  | 55kDa                                                     |
| <b>Uniprot ID</b>          | <a href="https://www.uniprot.org/entry/Q06413">Q06413</a> |
| <b>Dilution Range</b>      | WB: 1:500-2000, IHC-P: 1:50-200                           |
| <b>Expiration Date</b>     | 12 months from date of receipt.                           |



Western blot analysis of extracts of various cell lines.



Immunohistochemical analysis of formalin-fixed paraffin-embedded tissue.