



Product Datasheet

AQP0 antibody (orb135908)


Description	Rabbit polyclonal antibody to AQP0
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide of human MIP(NP_036196.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	28kDa
Uniprot ID	P30301
Dilution Range	WB: 1:500-1000, IHC-P: 1:100-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human st...



Immunohistochemical staining of rat brain...



Western blot analysis of extracts of mouse...