



Product Datasheet

Myosin IIa antibody (orb135951)


Description	Rabbit polyclonal to Myosin IIa
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Recombinant fusion protein of human MYH9(NP_002464.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	250kDa
Uniprot ID	P35579
Dilution Range	WB: 1:1000-2000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human st...



Western blot analysis of extracts of var...



Immunohistochemical staining of human es...