


## Product Datasheet

### MAGT1 Human Over-expression Lysate (orb1359750)

|                        |                                                                             |
|------------------------|-----------------------------------------------------------------------------|
| <b>Description</b>     | MAGT1 HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | MAGT1                                                                       |
| <b>Note</b>            | For research use only                                                       |
| <b>Uniprot ID</b>      | <b>Q9H0U3, A0A087WU53</b>                                                   |
| <b>Expiration Date</b> | 6 months from date of receipt.                                              |



Western  
blot  
validation  
of Left: Cell  
ly...