

Product Datasheet

NGFR antibody (orb135987)

Description	Rabbit polyclonal antibody to NGFR
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human p75NTR/NGFR(NP_002498.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45kDa
Uniprot ID	P08138
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.

