


Product Datasheet

PAK1 antibody (orb136039)

Description	Rabbit polyclonal antibody to PAK1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human PAK1(NP_002567.3).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	68KDa
Uniprot ID	Q13153
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, ICC: 1:50-200, IP: 1:20-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of mou...