



Product Datasheet

REG3A antibody (orb136041)


Description	Rabbit polyclonal antibody to REG3A
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human REG3A (NP_002571.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	15kDa
Uniprot ID	Q06141
Entrez	5068
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse stomach tissue.



Immunohistochemical staining of mouse brain tissue.



Immunofluorescence analysis of MCF-7 cells.