



Product Datasheet

PGC antibody (orb136069)


Description	Rabbit polyclonal antibody to PGC
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide of human II (PG II)(NP_001159896.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45-50kDa
Uniprot ID	P20142
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



IHC-P staining of human rectal tissue us...



Immunohistochemical analysis of formalin...



Western blot analysis of extracts of var...