

Product Datasheet

Substance P antibody (orb13610)

Description

Substance P antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human Substance P. Please contact us for the exact immunogen sequence. The peptide is available as orb374771.

Target

Substance P. Substance P antibody detects endogenous levels of total Substance

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

15 kDa

Uniprot ID

[P41539](#), [P06767](#), [P20366](#)

NCBI

[013998](#), [21](#), [054704](#), [11](#)

Entrez



[6863](#)

Dilution Range


IHC-P: 1:100-600, WB: 1:100-1000

Expiration Date


12 months from date of receipt.



IHC-P image of rat pancreas tissue using...



Western blot analysis of rat spinal cord...



Immunohistochemical staining of paraffin...