


Product Datasheet

MSH alpha antibody (orb136109)

Description	Rabbit polyclonal antibody to MSH alpha
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	A synthetic peptide of human POMC(NP_000930.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	29kDa
Uniprot ID	P01189
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...