



Product Datasheet

PPP2R4 antibody (orb136127)


Description	Rabbit polyclonal antibody to PPP2R4
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human Ser/thr-PP2A activator (PTPA)(NP_821067.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	31kDa/33kDa/36kDa/40kDa
Uniprot ID	Q15257
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...



IHC-P staining of rat brain tissue using...



Immunohistochemical analysis of formalin...