



Product Datasheet

PPP3R1 antibody (orb136130)

Description	Rabbit polyclonal to PPP3R1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human PPP3R1(NP_000936.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	15kDa
Uniprot ID	P63098
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100, IP: 1:20-50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemical staining of human pr...