

## Product Datasheet

### DPEP1 Human Over-expression Lysate (orb1361373)

|                        |                                                                             |
|------------------------|-----------------------------------------------------------------------------|
| <b>Description</b>     | DPEP1 HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | DPEP1                                                                       |
| <b>Note</b>            | For research use only                                                       |
| <b>Uniprot ID</b>      | <b>A0A140VJI3, P16444</b>                                                   |
| <b>Expiration Date</b> | 6 months from date of receipt.                                              |