



Product Datasheet

iNOS antibody (orb13614)


Description	iNOS antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human iNOS. Please contact us for the exact immunogen sequence. The peptide is available as orb374865.
Target	iNOS
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	131 kDa
Uniprot ID	P35228 , Q06518 , P29477
NCBI	31 , 41 , 000625 , 000616
Entrez	4843
Dilution Range	IHC-P: 1:100-600, WB: 1:100-1500
Expiration Date	12 months from date of receipt.



Western blot analysis of rat colon (lane...)



IHC-P staining of rat brain tissue using...



Immunohistochemical staining of guinea p...