



Product Datasheet

PRNP antibody (orb136166)


Description	Rabbit polyclonal antibody to PRNP
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human Prion Protein(NP_001073592.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20-30kDa
Uniprot ID	P04156
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of extracts of var...



IHC-P staining of mouse lung tissue usin...