



## Product Datasheet

### PTGIR antibody (orb136190)


|                            |                                                         |
|----------------------------|---------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal to PTGIR                              |
| <b>Reactivity</b>          | Human, Mouse, Rat                                       |
| <b>Conjugation</b>         | Unconjugated                                            |
| <b>Tested Applications</b> | IHC, WB                                                 |
| <b>Immunogen</b>           | Recombinant fusion protein of human PTGIR(NP_000951.1). |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide,50% glycerol,pH7.3.         |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.             |
| <b>Note</b>                | For research use only                                   |
| <b>Isotype</b>             | IgG                                                     |
| <b>Clonality</b>           | Polyclonal                                              |
| <b>MW</b>                  | 41kDa                                                   |
| <b>Uniprot ID</b>          | <b>P43119</b>                                           |
| <b>Dilution Range</b>      | WB: 1:500-2000, IHC-P: 1:50-200                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                         |



Immunohistochemical staining of human colon tissue showing PTGIR expression.



Immunohistochemical staining of human kidney tissue showing PTGIR expression.



Immunohistochemical staining of human stomach tissue showing PTGIR expression.