



Product Datasheet

RAB27A antibody (orb136214)


Description	Rabbit polyclonal antibody to RAB27A
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human RAB27A(NP_899057.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	27KDa
Uniprot ID	P51159
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:10-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human adrenal gland tissue.



Western blot analysis of extracts of various cell lines.



Immunohistochemical staining of human liver tissue.