



Product Datasheet

RNASE3 antibody (orb136247)

Description	Rabbit polyclonal to RNASE3
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human RNASE3(NP_002926.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	18kDa
Uniprot ID	P12724
Dilution Range	WB: 1:500-200, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human colon tissue.



Immunohistochemical staining of human stomach tissue.