



Product Datasheet

ompF antibody (orb13626)


Description	ompF antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human ompF. Please contact us for the exact immunogen sequence. The peptide is available as orb73071.
Target	ompF
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	42 kDa
NCBI	001005857 , 123441874
Dilution Range	P-ELISA: 1:500-1000, WB: 1:100-500, IHC-P: 1:100-500 (based on 0.5 mg/ml)
Expiration Date	12 months from date of receipt.



10% SDS-PAGE followed by a Western blot...



3&4, wild type(BW 25113) + kanamycin (0...



6&7, ΔOmpF (JW0912) + kanamycin (0xMIC,...