


## Product Datasheet

### p70S6Kb antibody (orb136260)

|                            |                                                                   |
|----------------------------|-------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal to p70S6Kb                                      |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                 |
| <b>Conjugation</b>         | Unconjugated                                                      |
| <b>Tested Applications</b> | IHC, WB                                                           |
| <b>Immunogen</b>           | Recombinant fusion protein of human p70 S6 Kinase 2(NP_003943.2). |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide,50% glycerol,pH7.3.                   |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.                       |
| <b>Note</b>                | For research use only                                             |
| <b>Isotype</b>             | IgG                                                               |
| <b>Clonality</b>           | Polyclonal                                                        |
| <b>MW</b>                  | 63kDa                                                             |
| <b>Uniprot ID</b>          | <b>Q9UBS0</b>                                                     |
| <b>Dilution Range</b>      | WB: 1:500-1000                                                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                   |



Western blot analysis of extracts of var...