


## Product Datasheet

# ARHGAP17 Human Over-expression Lysate (orb1362747)

|                        |                                                                                |
|------------------------|--------------------------------------------------------------------------------|
| <b>Description</b>     | ARHGAP17 HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | ARHGAP17                                                                       |
| <b>Note</b>            | For research use only                                                          |
| <b>Uniprot ID</b>      | <b>Q68EM7</b>                                                                  |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                 |



Western blot validation of Left: Cell ly...