



## Product Datasheet

### SDHA antibody (orb136291)


|                            |                                                        |
|----------------------------|--------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to SDHA                     |
| <b>Reactivity</b>          | Human, Mouse, Rat                                      |
| <b>Conjugation</b>         | Unconjugated                                           |
| <b>Tested Applications</b> | ICC, IF, IHC, WB                                       |
| <b>Immunogen</b>           | Recombinant fusion protein of human SDHA(NP_004159.2). |
| <b>Preservatives</b>       | PBS with 0.05% proclin300,50% glycerol,pH7.3.          |
| <b>Concentration</b>       | 1 mg/ml                                                |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.            |
| <b>Note</b>                | For research use only                                  |
| <b>Isotype</b>             | IgG                                                    |
| <b>Clonality</b>           | Polyclonal                                             |
| <b>MW</b>                  | 70KDa                                                  |
| <b>Uniprot ID</b>          | <b>P31040</b>                                          |
| <b>Dilution Range</b>      | WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:10-100      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                        |



Western blot analysis of extracts of var...



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...