


Product Datasheet

EAAT2 antibody (orb136308)

Description	Rabbit polyclonal to EAAT2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human EAAT2/EAAT2/SLC1A2(NP_004162.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65kDa
Uniprot ID	P43004
Dilution Range	WB: 1:500-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...