


## Product Datasheet

### Cyclin C (CCNC) Human Over-expression Lysate (orb1363095)

|                        |                                                                            |
|------------------------|----------------------------------------------------------------------------|
| <b>Description</b>     | CCNC HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | CCNC                                                                       |
| <b>Note</b>            | For research use only                                                      |
| <b>Uniprot ID</b>      | <b>Q7Z4L3</b> , <b>P24863</b>                                              |
| <b>Expiration Date</b> | 6 months from date of receipt.                                             |



Western blot validation of Left: Cell ly...