


## Product Datasheet

# FAM110C Human Over-expression Lysate (orb1364015)

|                        |                                                                               |
|------------------------|-------------------------------------------------------------------------------|
| <b>Description</b>     | FAM110C HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | FAM110C                                                                       |
| <b>Note</b>            | For research use only                                                         |
| <b>Uniprot ID</b>      | <b>Q1W6H9</b>                                                                 |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                |



Western blot validation of Left: Cell lysate