


## Product Datasheet

# TNFAIP8 Human Over-expression Lysate (orb1364022)

|                        |                                                                               |
|------------------------|-------------------------------------------------------------------------------|
| <b>Description</b>     | TNFAIP8 HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | TNFAIP8                                                                       |
| <b>Note</b>            | For research use only                                                         |
| <b>Uniprot ID</b>      | <b>O95379</b>                                                                 |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                |



Western  
blot  
validation  
of Left: Cell  
ly...