


## Product Datasheet

### LRRIQ4 Human Over-expression Lysate (orb1364438)

|                        |                                                                              |
|------------------------|------------------------------------------------------------------------------|
| <b>Description</b>     | LRRIQ4 HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | LRRIQ4                                                                       |
| <b>Note</b>            | For research use only                                                        |
| <b>Uniprot ID</b>      | <b>A6NIV6</b>                                                                |
| <b>Expiration Date</b> | 6 months from date of receipt.                                               |



Western  
blot  
validation  
of Left: Cell  
ly...