


## Product Datasheet

# LRRC4B Human Over-expression Lysate (orb1364439)

|                        |                                                                              |
|------------------------|------------------------------------------------------------------------------|
| <b>Description</b>     | LRRC4B HEK293T cell transient overexpression lysate (as WB positive control) |
| <b>Target</b>          | LRRC4B                                                                       |
| <b>Note</b>            | For research use only                                                        |
| <b>Uniprot ID</b>      | <b>A0A024R4I8, Q9NT99</b>                                                    |
| <b>Expiration Date</b> | 6 months from date of receipt.                                               |



Western  
blot  
validation  
of Left: Cell  
ly...