


Product Datasheet

PGC1 alpha antibody (orb13647)

Description	PGC1 alpha antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human PGC1 alpha. Please contact us for the exact immunogen sequence. The peptide is available as orb374773.
Target	PGC1 alpha
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	91 kDa
Uniprot ID	Q9UBK2 , Q9QYK2 , O70343
NCBI	013261 , 037393 , 31 , 11
Entrez	10891
Dilution Range	WB: 1:200-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of rat heart (lane...)