



Product Datasheet

BEST1 antibody (orb136488)


Description	Rabbit polyclonal to BEST1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human BEST1(NP_004174.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70kDa
Uniprot ID	076090
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human st...



Immunohistochemical staining of human li...



Immunohistochemical staining of rat brai...