


## Product Datasheet

### LCN1 Human Over-expression Lysate (orb1365287)

|                        |                                                                               |
|------------------------|-------------------------------------------------------------------------------|
| <b>Description</b>     | LCN1 HEK293T cell transient overexpression lysate<br>(as WB positive control) |
| <b>Target</b>          | LCN1                                                                          |
| <b>Note</b>            | For research use only                                                         |
| <b>Uniprot ID</b>      | <a href="#">A0A024R8D7</a> , <a href="#">P31025</a>                           |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                |



Western  
blot  
validation  
of Left: Cell  
ly...